

Table of contents

\\USER

<table border="1"><tr><td>Polimeni</td></tr></table>	Polimeni
Polimeni	

<table border="1"><tr><td><table border="1"><tr><td>Steffen</td></tr></table></td></tr></table>	<table border="1"><tr><td>Steffen</td></tr></table>	Steffen
<table border="1"><tr><td>Steffen</td></tr></table>	Steffen	
Steffen		

<table border="1"><tr><td><table border="1"><tr><td>260813-knee-openmsk</td></tr></table></td></tr></table>	<table border="1"><tr><td>260813-knee-openmsk</td></tr></table>	260813-knee-openmsk
<table border="1"><tr><td>260813-knee-openmsk</td></tr></table>	260813-knee-openmsk	
260813-knee-openmsk		

<table border="1"><tr><td><table border="1"><tr><td>wip_dess_CS_260902</td></tr></table></td></tr></table>	<table border="1"><tr><td>wip_dess_CS_260902</td></tr></table>	wip_dess_CS_260902
<table border="1"><tr><td>wip_dess_CS_260902</td></tr></table>	wip_dess_CS_260902	
wip_dess_CS_260902		

\\USER\Polimeni\Steffen\260813-knee-openmsk\wip_dess_CS_260902

TA: 2:34 min Coil Selection: Auto Voxel Size: 0.5×0.5×1.5 mm³ Acc:: 3.0 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	R71.1 A19.5 H89.4 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	80
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	160 mm
FOV Phase	100.0 %
Slice Thickness	1.50 mm
TR	15.19 ms
TE	4.90 ms
Averages	1
AutoAlign	---

Contrast - Common

TR	15.19 ms
TE	4.90 ms
Flip Angle	25 deg
Fat-Water Contrast	Water Excitation
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Multiple Series	Each Measurement

Resolution - Common

FOV Read	160 mm
FOV Phase	100.0 %
Slice Thickness	1.50 mm
Base Resolution	320

Resolution - Common

Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	CS
Total Factor	3.0
Reference Scans	Integrated
Acceleration Factor PE	2
Reference Lines PE	24
Acceleration Factor 3D	1
Reference Lines 3D	24
Reordering Shift 3D	0
Phase Partial Fourier	Off
Slice Partial Fourier	Off
Asymmetric Echo	Off
Elliptical Scanning	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	2D
Normalize	Off
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	R71.1 A19.5 H89.4 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	80
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	160 mm
FOV Phase	100.0 %
Slice Thickness	1.50 mm
TR	15.19 ms
Multi-Slice Mode	Sequential
Series	Interleaved

Geometry - AutoAlign

Slab Group	1
Position	R71.1 A19.5 H89.4 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	R71.1 A19.5 H89.4
R	71.1 mm

Geometry - AutoAlign

A	19.5 mm
H	89.4 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
CoilShim	Off
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	R71.1 A19.5 H89.4 mm
Orientation	Sagittal
Rotation	0.00 deg
A >> P	160 mm
F >> H	160 mm
R >> L	120 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	123.254122 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - PACE

Resp. Control	Off
---------------	-----

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Composing

Inline Composing	Off
------------------	-----

Inline - Open Recon

Algorithm	openmsk V0
Send original images	On
Compute thickness	On
Analysis	DESS
Segmentation model	DOSMA DESS
	bone/cartilage July 2024

Sequence - Part 1

Sequence Name	de
Dimension	3D
Excitation	Slab-sel.
RF Pulse Type	Normal
Gradient Mode	Fast
Flow Compensation	None
Bandwidth	260 Hz/Px
Asymmetric Echo	Off
Segments	1

Sequence - Part 2

Introduction	On
--------------	----

Sequence - Special

Interpolation Mode	Default
DL Recon	DL-CS-SR-CAIPI-P1
Run in Background	Off
Denoising Strength	0.0
CS Regularization	0.0010
CS Regularization 2	5.0
Gating States	0
Combining Mode	Addition
Scaling Factor for SE	1.0
Save Uncombined	On
Generate T2 Map	Off
Add Colorbar	Off

Sequence - Special

Use FA correction	Off
Symmetric SE	Off
PlayOutDiff	Off
Diff Mom. (R/P/S) 0	0.00 ms*mT/m
Diff Mom. (R/P/S) 1	0.00 ms*mT/m
Diff Mom. (R/P/S) 2	31.33 ms*mT/m

Sequence - Assistant

SAR Assistant	Off
---------------	-----